

Eight skincare packs (Tripaks) are included with your Omnilux clear-U device.

These skincare packs contain specially-formulated skincare to use in combination with Omnilux clear-U to maximize your results.



Follow these instructions on each light treatment day:

1. Cleanse your face using the Preparation Tripak.
2. Apply the Treatment Serum from the Treatment Tripak to your entire face.
3. Use Omnilux clear-U as directed for your 20-minute light treatment session. Do not remove Treatment Serum.
4. After the treatment, apply the Aftercare Tripak to your entire face.
5. Repeat above steps, alternating between blue and red light every 48 hours.

Please consult your User Guide for more information on your Omnilux clear-U.

Property of Photo Therapeutics, Inc.
Not to be duplicated or copied.

Product	What it does	Role in light therapy session
Preparation	Gently cleanses skin, removing makeup, dirt and debris. Improves circulatory flow and prepares and conditions the skin to accept the Treatment Serum.	The skin is better able to absorb light when it is clean and free from dirt and debris. When circulatory flow is increased from the Preparation combined with the Omnilux clear-U light, the active ingredients in the product are more readily absorbed.
Treatment	Rids the skin of debris and dirt that clog pores and reduces the apparent size of pores. Accelerates healing of irritated skin, encourages cell turnover and promotes overall healthier skin.	Works in combination with the Omnilux clear-U to improve nutrient exchange between cells. Uses an advanced transdermal delivery system to further improve the absorption rate of key ingredients that detoxify the skin. Healthy skin that is not irritated or inflamed is better able to absorb more light during the Omnilux clear-U treatment session.
Aftercare	Provides intensive hydration, improves elasticity and enhances the skin's rate of moisture retention. Repairs damage from the environment and improves the skin's immune system. Clarifies and evens skin tone.	Cells that are energized during the Omnilux clear-U session remain in a heightened state for up to 72 hours. Aftercare supplies continued nourishment to these cells in order to exploit the energized state of the cell.

Featured ingredients and benefits
<ul style="list-style-type: none"> • ROSEMARY LEAF EXTRACT and TEA TREE OIL, gentle astringents, deeply cleanse and unclog pores. • SALICYLIC ACID, APRICOT SEED AND APRICOT OIL, rich natural sources of AHA, gently exfoliate the skin. • ALOE BARBADENSIS LEAF EXTRACT and TOCOPHERYL ACETATE (Vitamin E) soothe the skin.
<ul style="list-style-type: none"> • Encapsulated PLANKTON EXTRACT promotes increased cellular metabolism. • WILLOW BARK EXTRACT unclogs pores, reduces the apparent size of pores and soothes skin by reducing areas of inflammation. • ZINC PCA repairs skin by helping to process the essential fatty acids that encourage healing. • RETINOL aids in collagen synthesis and promotes cellular renewal.
<ul style="list-style-type: none"> • SOY PROTEIN delivers deep hydration. • LECITHIN builds the lipid layer of cells and improves moisture retention. • ZINC PCA repairs skin by helping to process the essential fatty acids that encourage healing, and plankton extract improves the skin's immune system.