

Omnilux new-U Recommended Treatment Protocol

The new-U handheld LED device comes with 8 Tripaks of specially-formulated skincare that works in combination with the light to help improve periorbital wrinkles and promote healthy skin.

Follow these simple instructions for beautiful skin:

- Cleanse the face using the Preparation tripak
- Apply the Treatment serum from the Treatment tripak to your entire face
- Use Omnilux new-U as directed
- Apply the Aftercare tripak to the entire face over the Treatment serum.
- Wait at least 48 hours before your next new-U light treatment



We recommend using Photo Therapeutics' Firming Tripaks, as the efficacy of treatment has been proven in studies using this proprietary skincare product.

Is new-U safe and painless?

Yes. Subjects in our clinical trials have reported some skin redness after the first 1-2 treatments, but this is only temporary and will fade. No other adverse side effects were reported.



OMNILUX new-U™

for the effective treatment
of fine lines & wrinkles



PHOTO
THERAPEUTICS
world leader in light therapy

Contact info for USA:

2720 Loker Ave West Ste C, Carlsbad, California 92010
Tel: 760-607-0488 | Toll free: 800-743-8150 | Fax: 760-607-0288
E-mail: handheld@phototherapeutics.com

Contact info for ROW:

Unit 3 Bonehill Mews, Lichfield Street, Fazeley, Staffordshire B78 3QU
Sales: +44 (0) 1827 252 830 | Support: +44 (0) 161 925 5610
Facsimile: +44 (0) 1827 287 391
E-mail: salesinfo@phototherapeutics.com

www.phototherapeutics.com

The only clinically proven
hand-held device
FDA cleared for home use



Over 3300 dermatologists worldwide have been using Professional Omnilux products to successfully treat wrinkles and the effects of aging. Omnilux new-U™ uses the same LED technology but in a small and convenient home use package.

Omnilux new-U uses light technology in the form of light emitting diodes to:

- Stimulate natural anti-aging processes
- Improve skin tone and clarity

What is new-U?

Omnilux new-U™ is a handheld, Light Emitting Diode (LED) device for the treatment of periorbital wrinkles, more commonly known as crowsfeet. With the flip of a switch, its patented technology bathes the skin in alternating red and near infrared light. A study has also demonstrated improvements in skin tone and clarity. **Omnilux new-U is the first over-the-counter FDA cleared hand-held device for treatment of periorbital wrinkles.**

How does new-U work?

Omnilux combination phototherapy works with the body's own natural processes to counteract the effects of aging. The different colored lights target specific cells which are responsible for the synthesis and repair of the skin's supportive structures.

Who is suitable for new-U?

Those who have sun-damaged skin and/or fine lines and wrinkles can use new-U.



ACTUAL SIZE:
5"H x 2.5"W x 1.5"D, 2.9oz.



Photos courtesy Neil Sadick, New York NY

In an independent clinical study using Omnilux new-U, 7 out of 10 subjects reported a visible reduction in periorbital wrinkles and 8 out of 10 reported an improvement in skin clarity and smoothness after the recommended course of treatments.

What results can be expected from a series of treatments?

Best results will be seen between 4-8 weeks after treatment. The light stimulates natural processes that continue even after the treatment has stopped.

How many treatments will I need?

The course of treatments is 2x20 minute treatments per week for 4 weeks of alternate infrared and red light, but because the light will not damage your skin, there is no limit on the number of treatments you can have.

How long do the results last?

This depends on the age of your skin and the amount of photodamage. Skin rejuvenation is a dynamic process, and natural expression lines (i.e. lines caused by laughing, smiling, frowning and crying) will gradually return in the course of time. Since we can't stop the aging process; we recommend follow-up treatments (typically 2 to 4 treatments) every month dependent upon your lifestyle.

Is new-U all I need for my skin?

Omnilux new-U is not a substitute for total skincare. We suggest you speak with your skincare specialist to find the best combination for your skin.

order today...

USA orders: Call 800-743-8150 or email handheld@phototherapeutics.com

ROW orders: Call +44 (0) 1827 252 830 salesinfo@phototherapeutics.com