

Eight skincare packs (Tripaks) are included with your Omnilux new-U device.

These skincare packs contain specially-formulated skincare to use in combination with Omnilux new-U to maximize your results.



Follow these instructions on each light treatment day:

1. Cleanse your face using the Preparation Tripak.
2. Apply the Treatment Serum from the Treatment Tripak to your entire face.
3. Use Omnilux new-U as directed. Do not remove Treatment Serum.
4. After the treatment, apply the Aftercare Tripak to your entire face.
5. Repeat above steps, alternating between red and infrared light every 48 hours.

Please consult your User Guide for more information on your Omnilux new-U.

Property of Photo Therapeutics, Inc.
Not to be duplicated or copied.

Product	What it does	Role in light therapy session
Preparation	Cleanses skin, removing makeup, dirt and debris. Improves circulatory flow and prepares and conditions the skin to accept the Treatment Serum.	The skin is better able to absorb light when it is clean and free from dirt and debris. When circulatory flow is increased from the Session Preparation combined with the Omnilux new-U light, the active ingredients in the Treatment product are more readily absorbed.
Treatment	Increases rate of cellular energy and supplies nutrients that are required to produce elastin to cells that are energized during the session. Provides intensive hydration and improves moisture retention. Promotes immediate lifting and improves elasticity. Removes impurities and detoxifies the skin.	Works in combination with the Omnilux new-U to improve metabolism. Feeds energized cells the nutrients required to produce collagen. Uses an advanced transdermal delivery system to further improve the absorption rate of key ingredients. Oxygen free radicals are a natural by-product of elastogenesis, the process of manufacturing elastin. Antioxidants in the Session Treatment product work to reduce oxidative stressors and rid the skin of free radicals without disrupting the natural chain of events.
Aftercare	Provides continued nourishment to the skin by supplying the nutrients used in producing collagen and elastin. Enhances oxidative stress management and improves the skin's photo-immune system, enabling better defense against the environment.	Cells that are energized during the Omnilux new-U session remain in a heightened state for up to 72 hours after the session. Aftercare supplies continued nourishment to these cells in order to exploit the energized state of the cell and encourage maximum production of collagen and elastin. The body produces collagen and elastin through an oxidative process from which free radical oxygen molecules are a by-product. The fragile network of collagen and elastin fibres that is developing requires enhanced protection from these free radicals, as well as other environmental insults. By preserving collagen and elastin from previous Omnilux new-U sessions, results are able to accumulate and build over multiple sessions.

Featured ingredients and benefits

- ROSEMARY LEAF EXTRACT and TEA TREE OIL, gentle astringents, deeply cleanse and unclog pores.
 - APRICOT SEED AND APRICOT OIL, both rich natural sources of AHA, gently exfoliate the skin.
 - ALOE BARBADENSIS LEAF EXTRACT and TOCOPHERYL ACETATE (Vitamin E) soothe the skin.
-
- COQ10 enhances cellular metabolism.
 - Encapsulated SOLUBLE RETINOL increases the rate of collagen and elastin synthesis.
 - SWEET ALMOND and ANISE EXTRACTS form a tensor film that provides an immediate lift and prevents expression lines caused by muscle contractions.
 - CENTELLA ASIATICA and CORNFLOWER EXTRACTS naturally increase skin thickness and improve firmness.
 - PASSION FLOWER, BLACKCURRANT and GRAPE LEAF EXTRACTS dramatically improve elasticity and inhibit the production of enzymes that break down elastin.
 - CARNITINE combines with ALGAE EXTRACT to act as a powerful oxidative reducing agent that battles free radicals.
-
- SODIUM PCA, MAGNESIUM PCA and ZINC PCA nourish the skin with vital precursors to collagen and elastin.
 - Photo-Activated VANILLA TAHITENSIS FRUIT EXTRACT and PLANKTON EXTRACT with LECITHIN leverage a liposomal delivery system to repair oxidative stress and DNA damage from sun exposure.
 - ESTER OF VITAMIN C acts as a powerful antioxidant and inhibits the appearance of sun spots and hyperpigmentation.
 - TOCOPHERYL ACETATE (Vitamin E) and ESTER OF VITAMIN C act as powerful antioxidants, ridding the skin of free radicals.